



www.sibelmed.com

DISPOSABLE OTOSCOPE SPECULUM

MODEL	HG4	HG2	HM4	HM2	R4	R2	R4L3	R2L3
DIAMETER	4 mm	2.5 mm	4 mm	2.5 mm	4 mm	2.5 mm	4 mm	2.5 mm
COMPATIBLE	<p>Heine (Diagnóstico, K100, Beta 100)</p> <p>Kawe (Eurolight)</p>	<p>Heine (Diagnóstico, K100, Beta 100)</p> <p>Kawe (Eurolight)</p>	<p>Heine (Kappa, Beta, Alpha, MiniLux, Mini FibraLux, Mini 2000, K180, Beta 200, Mini 3000)</p> <p>Kawe (piccolight)</p> <p>Riester (Ri-Mini, Pen-Scope)</p> <p>Welch Allyn (serie 524)</p>	<p>Heine (Kappa, Beta, Alpha, MiniLux, Mini FibraLux, Mini 2000, K180, Beta 200, Mini 3000)</p> <p>Kawe (piccolight)</p> <p>Riester (Ri-Mini, Pen-Scope)</p> <p>Welch Allyn (serie 524)</p>	<p>Riester (Ri-Scope L1 y L2, Ri-Star)</p>	<p>Riester (Ri-Scope L1 y L2, Ri-Star)</p>	<p>Riester (Ri-Scope L3)</p>	<p>Riester (Ri-Scope L3)</p>

Available in 2 formats (from 1,000u.)



Dispenser Box of 250u.



20 Strips of 50u.

Use products certificated in accordance with 93/42/EEC Directive in order to guarantee the patient security as well as the quality of the test.

SIBEL, S.A

Rosselló, 500 08026 Barcelona (España)

T.: +34 93 436 00 07 Fax: +34 93 436 16 11

export@sibelmed.com

SIBEL, S.A. belongs to

Sibel GROUP